## UHT Product Specification for Lakeland UHT Semi Skimmed Stick 240 x 10ml Status: Live Revision: 9

General Information					
Product	Lakeland UHT Semi Skimmed Stick 240 x 10ml				
Resource Number	0874721				
Production Site	Killeshandra				
Inner Barcode	None	Outer Barcode	05099869747213		
Gross Weight (Kg)	2.711	Net Weight (Kg)	2.477		
Palletisation	Cases per layer: 30 Layers on pallet: 8 Qty on pallet: 240	Outer Case Dimensions	Ext 193x193x135		
Created	17/08/2016	Created By	Louise McKiernan		
Last Updated	16/11/2023 11:05:58	Last Updated By	Nuala McBride		

Description	UHT Semi Skimmed Milk Portions
Catalogue Name	
Use/Benefits	Ideal for use in tea or coffee, Ambient storage, Portion control.
Packaging	240 x 10ml aseptic stick
Shelf Life	7 Months
Ingredients	Semi Skimmed Milk
Allergy Advice	For allergens, see ingredients in <b>bold</b>
Allergy Advice (only for	
Custom)	

Nutrient  nergy at of which saturates arbohydrate of which sugars	per 100 ml 201kJ / 48kcal 1.5g 1g		
nergy at of which saturates arbohydrate	201kJ / 48kcal 1.5g		
at of which saturates arbohydrate	1.5g		
of which saturates arbohydrate			
arbohydrate	1g		
f which gugara	4.9g		
	4.9g		
rotein	3.6g		
alt	0.1g		
It content is exclusively	vidue to the presence of natu	rally occurring sodium	
Once opened use immediately			
Store in a cool dry place below +20°C. Once opened use immediately			
our Smooth/ White			
Clean Taste			
Absent			
is			
1.6%	1.6%		
Sterility is assured	Sterility is assured by incubation of representative batch samples at 32°C for 3 days prior to analys enzymatic bioluminescence.		
		nder aseptic conditions.	
-	case only)		
DD/MM/YY			
9:00	es our Smooth/ White Clean Taste Absent  Sterility is assured enzymatic biolum Sterilised by Ultra Y DDD (outercase)	es our Smooth/ White Clean Taste Absent  1.6%  The product is sterilised by UHT treatment and Sterility is assured by incubation of representation enzymatic bioluminescence.  Sterilised by Ultra Heat Treatment. Packed un Y DDD (outercase only)	pes our Smooth/ White Clean Taste Absent  1.6%  The product is sterilised by UHT treatment and aseptically packed. As such it is con Sterility is assured by incubation of representative batch samples at 32°C for 3 days enzymatic bioluminescence.  Sterilised by Ultra Heat Treatment. Packed under aseptic conditions. Y DDD (outercase only)

Legislation	The product meets all relevant UK & EC legislation
Other Information	
To view a previous vers	sion of this document click on the links below.