

## Tork Blue Couch Roll Advanced



Article	152250
Roll Length	56.1 m
System	C1 - Couch roll system
Roll Width	48 cm
Roll Diameter	12 cm
Number of Sheets	165
Sheet Length	34 cm
Core Inner Diameter	4.4 cm
Ply	2
Print	No
Embossing	Yes
Colour	Blue

Protect your patient couch with the disposable Tork Couch Roll Advanced. The perforated couch roll is ideal to maintain a high level of hygiene in the healthcare environment. Rolls can be used with the Tork Couch Roll Dispenser.

[www.tork.co.uk](http://www.tork.co.uk)

Better resistance and higher comfort due to the double ply

Perforation for easy tear-off and consumption control

## Shipping data

	Consumer unit	Transport unit	Pallet
<b>EAN</b>	7310791230375	7310791230320	7322540716146
<b>Pieces</b>	1	9	162
<b>Consumer units</b>	-	9	162
<b>Height</b>	480 mm	480 mm	1121 mm
<b>Width</b>	120 mm	360 mm	1000 mm
<b>Length</b>	120 mm	360 mm	1200 mm
<b>Volume</b>	6.9 dm <sup>3</sup>	62.2 dm <sup>3</sup>	1.1 m <sup>3</sup>
<b>Net weight</b>	969 g	8.7 kg	157 kg
<b>Gross weight</b>	1001 g	9.1 kg	163.8 kg
<b>Packing material</b>	Plastic	Plastic	-

## Choose a dispenser



129183

129184

## Alternative products



152225

150250

150225

## Product certificates



## Contact

Essity UK Ltd

Southfields Road  
Dunstable  
Bedfordshire  
LU6 3EJ