

Blue Absorbent Mat 500 x 580mm
Product Specification

- Product:** Blue Absorbent Mat 500 x 1000mm
- Description:** Blue Absorbent Mats are an absorbent mat used underfoot in operating theatres to absorb liquid spillage. It has a non-woven top and a blue high grip perforated PE underside. The perforations allow water to be drawn up from under the mat thus ensuring good grip. The mat is 580mm long and 1000mm wide. It is made using custom raw materials.
- Structure:** Blue Absorbent Mat 500 x 1000mm Mats consist of:



- Absorbency:** The standard absorbency for Sirane Dri- Fresh[®] 3000 PF Blue Theatre Mats is 3000cc/m²
 Nominal absorbency per pad: 1740cc/ linear meter
 Minimum absorbency tolerance is 25% below the nominal.
**Due to processing of material and lamination properties the absorbency of Sirane products will on occasion be below the nominal figure.
 Testing should be conducted to ensure that the pad and required absorbency is suitable for the required application.*

- Dimensions:** 580mm x 1000mm.
- Width Tolerance:** ± 5mm
- Length Tolerance:** ± 30mm
- Weight:** Nominal individual pad weight: 110 g
- Thickness:** Sirane do not specify thickness, only absorbency.
- Packaging:** Boxed in blue, food-grade liners.
- Palletising:** Boxes shrink-wrapped on standard 1200mm x 1000mm pallet.
- Labelling:** Individual box label includes product description, dimensions, colour, traceability code, customer specified product code, order number, box quantity and box number. Pallet label includes customer name, delivery address, product description, dimensions, colour, traceability code, customer specified product code, box quantity, pallet quantity and shipping date.

Sirane Dri-Fresh[®] 3000 PF Blue Theatre mats are manufactured in BRCGS Packaging Materials, ISO13485 and ISO 9001:2015 accredited facilities at the Sirane Development Centre in Telford.

Certification for food contact, ISO and BRCGS accreditation are available on request.

Disclaimer: This specification is intended as a general guide only. It is recommended that full trials are carried out to ensure the suitability of Sirane Dri-Fresh[®] 3000 PF Blue Theatre mats for food packaging applications.